

MATERIAL SAFETY DATA SHEET

SECTION I: IDENTIFICATION

<u>Catalog Number(s)</u>	<u>Chemical/Common Name</u>	<u>Formula</u>	<u>Hazardous</u>
3000-28001	BiokitHSV-2 Rapid Test 20 Test		No

Biokit USA, Inc.
101 Hartwell Avenue
Lexington, MA 02421

MSDS Coordinator
Telephone Number: (781) 861-4066

Date Prepared: 2/16/06

Emergency Phone: (800) 926-3353

SECTION II: HAZARDOUS INGREDIENTS

Pursuant to OSHA's *Hazardous Communications Standard* 29 CFR 1910.1200, these components are not considered to contain hazardous materials. They do, however, contain 0.1% sodium azide as a preservative for which we are providing a Material Safety Data Sheet.

Sodium Azide **CAS-Number:** 26628-22-8

<u>Trade Names/Synonyms</u>	<u>Chemical Family</u>	<u>Molecular Formula</u>	<u>Molecular Weight</u>
Azlum; Azide; Kazoe	Inorganic Salt	NaN ₃	65.01

SECTION III: PHYSICAL DATA

<u>Appearance</u>	<u>Melting Point</u>	<u>Specific Gravity</u>	<u>Solubility In Water</u>
Colorless Crystalline Solid	275°C (527°F) decomposes	1.8 (water = 1.0)	400mg/mL

SECCION IV: FIRE AND EXPLOSION HAZARD DATA

Extinguishing Media

Sand, sodium chloride, sodium carbonate, dry chemical powder; do not use water.

Special Firefighting Procedures

Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

Unusual Fire and Explosion Hazards

Container explosion may occur under fire conditions. Emits toxic fumes under fire conditions.

SECTION V: REACTIVITY DATA

Stability

Unstable

Conditions to Avoid

- Avoid contact with metals and acids
- Explodes when heated
- May be shock sensitive

Incompatibilities

- Acid chlorides
- Nitrogen oxides
- Azide reacts with many heavy metals such as lead, copper, mercury, silver and gold to form explosive compounds. Copper and lead azides are more sensitive than nitroglycerin. Azide reacts with metal halides to give a range of metal azide halides, many of which are explosive. Incompatible with chromyl chloride, hydrazine, bromine, carbon disulfide, dimethyl sulfate, dibromomalonitrile.

Hazardous Combustion or Decomposition Products

Toxic fumes of Nitrogen oxides

Hazardous Polymerization

Will not occur

SECTION VI: HEALTH HAZARDS

BK HSV-2 is intended for IN VITRO diagnostic use only.

Carcinogen Status

None....Sodium azide is a strong skin and mucous membrane irritant and poison. Contact may produce skin burns. Poisoning from inhalation or ingestion may result in anesthetic symptoms such as headache, dizziness, faintness, somnolence, and unconsciousness. Serious poisonings may be fatal.

First Aid

- In case of exposure, obtain medical attention immediately. If ingested, give large quantities of water and induce vomiting.
- In case of skin contact, flush with copious amounts of water for at least 15 minutes. Remove contaminated clothing and shoes and call a physician. If inhaled, remove to fresh air. If breathing becomes difficult, call a physician.
- In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating the eyelids with fingers. Call a physician.

SECTION VII: PRECAUTIONS FOR SAFE HANDLING AND USE

Ventilation

Provide local exhaust ventilation system to meet permissible exposure limits.

Respiration

High Levels...Supplied-air respirator with a full facepiece, helmet or hood. Self-contained breathing apparatus a full facepiece.

Firefighting

Self-contained breathing apparatus with full facepiece operated in a pressure demand or other positive-pressure mode.

Clothing

Protective clothing not required. Avoid repeated or prolonged contact with this substance.

Gloves

Employee must wear appropriate protective gloves to prevent contact with this substance.

Eye Protection

Employee must wear splash-proof or dust-resistant safety goggles and a faceshield to prevent contact with this substance.

SECTION VIII: CONTROL MEASURES

Contain spill with absorbent material and place in appropriate container. Wash and disinfect spill site after clean-up.

SECTION IX: ADDITIONAL COMMENTS

Store 2°-8°C

Employers should use this information as a supplement to other information gathered by them and should make independent judgment of the suitability of this information to ensure proper use and to protect the health and safety of employees. This information is provided without warranty and any use of this product not in conformance with this Material Safety Data Sheet, or in combination with any other product or process is the responsibility of the user. Biokit USA, Inc. shall not be held liable for any damage resulting from the handling or use of this product.

MANUAL - SAFETY
DATE REVIEWED _____
INITIALS _____

LOCATION/ DEPARTMENT

FORM NO. _____

ISSUE DATE _____
SUPERSEDES MSDS
FORM NO. _____ DATED _____